

GLUCOMETER CHECKOFF

COMPLETED

Applies appropriate PPE	
Gathers equipment/supplies	
Informs patient of procedure	
Cleanses testing site (distal fingertip)	
Turns on device by inserting testing strip	
Applies lancet to patients fingertip	
Allows testing strip to draw in appropriate amount of blood	
Places used lancet in sharps container	
Wipes testing site clean and holds or instructs patient to hold direct pressure	
Applies bandaid to patients fingertip to stop any bleeding	
Disposes of used testing strip	
Records results of glucometer testing in emscharts	
Cleans glucometer per manufacturer recommendations	
Places glucometer back in case	
Demonstrates glucometer checkoff per the truck check off form	

By signing this both the employee and the skills evaluator verify competency with glucometer

Name of employee_____

Signature of employee_____

Date Signed_____

Name of skills evaluator_____

Signature of skills evaluator_____

Date Signed_____