



STAIR CHAIR USAGE

Name: _____ CERT LEVEL: _____

Date: _____

COMPLETED	SKILL
P / F	Obtains stair chair and opens the chair
P / F	Ensures the stair chair has locked into position
P / F	Checks the stair chair to ensure straps are in working order and stair chair is clean and in good working order
P / F	Explains the procedure to the patient and answers any questions appropriately
P / F	Assist the patient to the stair chair and positions them properly in the chair
P / F	Secures the patient correctly with all straps being secured
P / F	With a partner approaches the stairs with the patients back facing the stairs
P / F	Checks with partner to ensure they are ready to initiate the lift
P / F	On a count of 3 initiates the lift by leaning back the chair and pulling the chair towards them
P / F	Ensures the chair maintains a centered position on the stairs with each lift
P / F	Once at the top of the stairs, ensures the chair is on level ground away from the stairs with wheels locked
P / F	Assist the patient off of the stair chair and to a standing position. Properly supports the patient to prevent fall risk
P / F	Gathers and puts away all equipment in designated area

Comments:

EVALUATOR NAME:

SIGNATURE: